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- 2) Azuma H, Hirayama J, Akino M, Miura R, Kiyama Y, Imai K, Kasai M, Koizumi K, Kakinoki Y, Makiguchi Y, Kubo K, Atsuta Y, Fujihara M, Homma C, Yamamoto S, Kato T, Ikeda H. Reduction in adverse reactions to platelets by the removal of plasma supernatant and resuspension in a new additive solution (M-sol). *Transfusion* 49: 214-218, 2009
- 3) Yanagisawa R, Shimodaira S, Kojima S, Nakasone N, Ishikawa S, Momose K, Honda T, Yoshikawa K, Saito S, Tanaka M, Nakazawa Y, Sakashita K, Shiohara M, Akino M, Hirayama J, Azuma H, Koike K. Replaced platelet concentrates containing a new additive solution, M-sol: safety and efficacy for pediatric patients. *Transfusion* 53: 2053-2060, 2013
- 4) Berseus O, Boman K, Nessen SC, Westerberg LA Risk of hemolysis due to anti-A and anti-B caused by the transfusion of blood or blood components containing ABO-incompatible plasma. *Transfusion*, 53:114S-123S, 2013
- 5) 佐々木 大、小砂子 智、小宮山祥光、鈴木 光、浦野慎一、中野月子、清水哲夫、神谷 忠、平沼 隆明、西岡 克郎、伊田八洲雄. 血小板の洗浄・保存液の比較検討. *日本輸血学会誌* 47: 777-782, 2001
- 6) Hirayama J, Azuma H, Fujihara M, Homma C, Yamamoto S, Ikeda H. Storage of platelets in a novel additive solution (M-sol), which is prepared by mixing solutions approved for clinical use that are not especially for platelet storage. *Transfusion*, 47: 960-965, 2007
- 7) 平山順一、東 寛、藤原満博、秋野光明、本間稚広、山本定光、加藤俊明、池田久實. 市販輸液の混合物である新たな洗浄置換液 (M-sol) による血小板の保存. *日本輸血細胞治療学会誌* 54: 17-22, 2008
- 8) Oikawa S, Sasaki D, Kikuchi M, Sawamura Y, Itoh T. Comparative in vitro evaluation of apheresis platelets stored with less 100% plasma versus bicarbonated Ringer's solution with less than 5% plasma. *Transfusion*, 53:655-660, 2013

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